

www.fluidinova.com

Nanocrystalline apatite: from fundamentals to bone substitute materials Industrial and business point of view



Paulo Quadros, PhD

Head of Product & Business Development

NEWGEN Workshop, Nantes, May 2014

About us



🛯 米 🏊 🛯 🖌 🖛 🚃 🖉





Hydroxyapatite Nanoparticles

Fluidinova is a specialized manufacturer of synthetic nanocrystalline hydroxyapatite material, which is commercialized under the brand name nanoXIM®. Using an innovative proprietary technology – NETmix® – Fluidinova produces a highly pure, single phase nanoHydroxyapatite. nanoXIM® is supplied worldwide in different forms to manufacture Medical Devices, for Personal Care Product formulations, R&D activities and many other applications.

Follow us on LinkedIn



NEWS - Check out the new video of nanoXIM - The Enamel Repair Ingredient.

Download product brochure

Business areas





for medical devices



for personal care



High-quality hydroxyapatite supplied to medical devices manufacturers.

Several clients using nanoXIM•HAp for manufacturing bone substitutes with CE mark in Europe and FDA approval in the USA.

Bone substitutes are the main nanoXIM Medical market, with special emphasis on the injectable bone substitutes devices.

nanoXIM possesses excellent properties of bone regeneration and biocompatibility.

Personal Care segment of the nanoXIM product

High-quality hydroxyapatite supplied for cosmetic manufacturers, with special focus on dentifrices (toothpastes, mouthwashes, whitening products...).

nanoXIM nanohydroxyapatite presents good enamel remineralization potential and high occlusion efficacy of dentine tubules, making it a solution for teeth hypersensitivity.

Tooth surface smoothness and additional whitening are a plus provided by nanoXIM.



for several applications

R&D applications of the nanoXIM product in a broad applications fields

High-quality hydroxyapatite material for R&D projects like bio-separation, pharma, drug delivery and polymers functionalisation.





nanoHAP industrial manufacturing

















R&D projects:

- NEWGEN: New Generation Biomimetic and Customized Implants for Bone Engineering, COST EU RTD (2013-2017)
- IMCOSS: Injectable Medical Ceramics for Bone Repair and Augmentation, FP7 EU (2012-2014)
- NanoForBone: Development and scale-up of 3D structures and coatings for medical applications based on nanosized hydroxyapatite, Adl Portugal (2009-2012)





Applications, Products & Studies





NANOXIM Medical

APPLICATIONS >

nanoXIM•Medical materials are precursors of high quality medical devices biomaterials used in:





Trauma & Orthopedics

Filling bone voids or gaps not intrinsic to the stability of the bony structure caused by trauma, disease or related surgical procedures reducing significantly the use of autologous grafts.

nanoXIM•Medical pastes are suitable to be injectable, moldable and easily adaptable to bone defect shape, allowing a high bone-implant interface.



Periodontal & Maxillofacial

Periodontal defects filling, alveolar bone filling, sinus base augmentation and elevation, reconstruction of mandibular cyst defects and voids after tooth socket extractions.



Otology

Components in otology implants and cements suitable for middle ear surgeries. Innovative otology products enable simpler surgical procedures, reducing recovery times, thus improving patient's quality of life.



Pharma

The high surface area and biocompatibility of nanoXIM spray-dried powders makes them particularly suitable for the development of new drug delivery systems.

< PRODUCTS



nanoXIM•Medical Pastes

> nanoXIM+HAp100 series are water based nano-hydroxyapatite pastes specially recommended for medical devices manufacturing such as bone injectable substitutes for bone repair and reconstruction.



ADVANTAGES

- > Highly osteostimulative
- > Highly resorbable material replaced by new bone during the healing process
- > Optimal defect filling due to pasty consistency
- > 100% synthetic and safe material
- > High surface area





Product Reference	Description
nanoXIM+HAp102	15.0 ±1.0% wt. hydroxyapatite nanoparticles aqueous paste
nanoXIM•HAp103	30.0 ±3.0% wt. hydroxyapatite nanoparticles aqueous paste

nanoXIM•Medical Powders

> nanoXIM+HAp200 series are spray-dried hydroxyapatite powders used as precursors of porous granules and blocks scaffolds for bone repair and reconstruction.



ADVANTAGES

- > 100% synthetic and safe material
- > High surface area and porosity
- > Nanostructured micron sized material
- > High biocompatibility
- > Narrow particle size distribution





Product Reference	Description
nanoXIM•HAp202	5.0 ± 1.0 μm d ₅₀ hydroxyapatite spray-dried powder
nanoXIM•HAp203	10.0 ± 2.0 µm d ₅₀ hydroxyapatite spray-dried powder





Studies



In vitro study with nanoXIM paste



HMSCs cultured in the absence (control) and in the presence of 10 μg/mL of nanoXIM•HAp102 at days 7, 14 and 21. (A) Histochemical staining for ALP and (B) collagen. The synthesis of ALP is observed by the presence of a brown to black stain and collagen by a pink to red stain. (C) Confocal laser scanning microscopy of HMSCs cultured for 14 days in the absence (control) and in the presence of 10 μg/mL of nanoXIM•HAp102. F-actin cytoskeleton organization is stained in green and nuclei in red.

• Human mesenchymal stem cells (HMSCs) are considered an excellent candidate for cell-based therapies since they are able to differentiate in a wide variety of cell types, including the osteogenic lineage;

• HMSCs exposed to 10 µg/mL of nanoXIM-HAp102 showed increased alkaline phosphatase (ALP) (A) and collagen synthesis (B) and a typical cytoskeleton organization (C);

• These results demonstrate that nanoXIM-HAp102 formulation is able to modulate the HMSCs osteoblastic differentiation, being a promising product to be used in bone regeneration applications.

Source: E. Pires, "Effect of the nanohydroxyapatite Formulation NanoXIM.HAp102 on the Proliferation and Osteogenic Differentiation of Human Bone Mesenchymal Stem Cells", Integrated MSc Thesis in Bioengineering, Faculty of Engineering, University of Porto (2013).



In vitro study with nanoXIM powder



MG63 osteoblast-like cells cultured on micro and nanoscale HA substrates. (A) Cell proliferation for 3 and 6 days estimated by resazurin assay; (B) Confocal microscopy images of cells cultured for 6 days on the two different materials. Data kindly provided by Marta Laranjeira.

• The substrates produced using nanoXIM-HAp202 showed **improved cell viability and proliferation** when compared with microscale hydroxyapatite substrates (A);

• Both substrates present well spread cells with typical morphology. However, nanoscale substrates possess a **higher number of cells** adhered on the material (B);

• Nanostructured hydroxyapatite substrates produced with nanoXIM-HAp202 provided a more adequate environment for bone regeneration.

Source: Laranjeira MS, Fernandes MH, Monteiro FJ. Innovative macroporous granules of nanostructured-hydroxyapatite agglomerates: Bioactivity and osteoblast-like cell behaviour. Journal of Biomedical Materials Research Part A. 2010;95A(3):891-900.





Next Generation Products



nanoXIM blocks and granules





nanoXIM blocks and granules



Granules 1-2 mm





FLUIDINOVA, SA

Tecmaia - Rua Eng^o Frederico Ulrich, n^o 2650 4470-605 Moreira da Maia PORTUGAL 20 119 746

www.fluidinova.com

nanoxim@fluidinova.com

