MP1301- NEWGEN WG1

Nanocrystalline apatite: from fundamentals to bone substitute materials

Christèle Combes (INPT, CIRIMAT, France) Margarida Almeida (University AVEIRO, CICECO, Portugal) Paulo Quadros (Fluidinova, Portugal)











Nanocrystalline apatites:

from fundamentals to bone substitute materials

<u>Christèle Combes</u>, Christophe Drouet, Sophie Cazalbou, Stéphanie Sarda, Christian Rey

University of Toulouse, CIRIMAT, UPS-INPT-CNRS, ENSIACET, Toulouse, France

christele.combes@ensiacet.fr



NEWGEN 1st Workshop, 7th-8th may 2014, School of Dentistry, Nantes, France



Biological mineralisation properties

- Mineral elements and structures quite similar Silicates Calcium carbonates Calcium phosphates
- <u>Adaptable and elaborate architecture</u> (growth and autorepair)
- Resistance and lightness
- Saving energy (low temperatures, recycling, renewal)

A model for material scientists !



Characteristics of biological mineralisations (vertebrates)

Bone

Tooth (enamel)

≠ Architecture

Composition of tissue

• 50-85 % mineral

• 97-99 % mineral

Crystal characteristics

- Nanocrystals

 (20-60 x 10-20 x 2-5 nm)
 Platelets
 High specific surface area
- large crystals

Needles low specific surface area



Introduction

16 Å

Structure and composition of bone tissue/mineral







Synthesis, properties and applications of nanocrystalline apatites analogous to bone mineral



Nanocrystalline apatites analogous to bone mineral <u>Synthesis</u>

 $Ca_{10-x} \square_{x} (PO_{4})_{6-y} (HPO_{4}, CO_{3})_{y} (OH)_{2-z} \square_{z}$

 $2x = y+z; x \le 2$

Ca $(NO_3)_2$, 4 H₂O $(NH_4)_2$ HPO₄ ambient T physiological pH



 \approx 11x19 nm (1d maturation)

Nanocrystalline apatites properties:

- Maturation (h, d, m)
- Ion exchanges
- Protein adsorption



X-ray diffraction characterisation



Estimated crystallite size for natural and synthetic apatites

[Eichert, PhD thesis, INP Toulouse 2001]

Sample	L(002) ± 3 Å (length)	L(310) ± 3 Å (width/thickness)
Chicken bone	207	66
Rabbit bone	190	Not evaluated
Adult human cortical bone	213	68
Synthetic apatite matured for 3 months	282	72



FTIR spectroscopy analysis





Apatites maturation

HPO₄²⁻ = important marker of apatite maturation



Synthetic apatites

Biological apatites (chicken bone)



HPO₄²⁻ **** with maturation/ageing

[Cazalbou *et al., J. Mater. Chem.,* 2004] [Cazalbou *et al., J. Bone Miner. Metab.,* 2004]



Ion exchange properties of nanocrystalline apatites

Cation exchanges Sr²⁺/Ca²⁺



fast and reversible exchanges exchange ability **\u00e4** with maturation

Anion exchanges CO₃²⁻/HPO₄⁻



Surface divalent ions involved in these ion exchanges

[Cazalbou et al., J. Mater. Chem., 2004]



Spectroscopic characterisation

Role of water ?

FTIR spectroscopy

³¹P NMR



Wet samples exhibit

fine structural details which disappear on drying



Drying:

- Loss of water
 - ⇒ disordered ionic layer
 - ⇒ alteration of short range order
 - Non-apatitic environments

Very fragile surface hydrated layer structure

Some analogy but **≠ OCP and DCPD**

[Eichert et al., Bioceramics 16, 2004]



Nanocrystalline apatite properties

Adsorption properties of nanocrystalline apatites



Adsorption properties vary according to the surface layer characteristics

■ ● adsorption
□ ○ desorption

B: stoechiometric HAP



Specific structural characteristics of apatite nanocrystals

Structured hydrated layer at the surface of nanocrystals





Consequences for biomaterials ?





Nanocrystalline apatites analogous to bone mineral



Biocompatibility, biomimetism High surface reactivity

Bone tissue engineering

- Low temperature bioceramics (biomimetic cements, SPS, gel casting, ...

Bioactivity Resorbability



Associations with drugs (antibiotic, anticancer, growth factor, biologically active ions ...) for sustained release system

- Nanocomposites (nanocrystalline apatites + resorbable polymers)



Ceramic prepared by SPS (Cryo SEM-FEG)



tem cells on nanostructured coating

(ceramics or metal)



+ association with biologically active molecules *Histological cross section : neoformed bone (dark blue)*

Istological cross section : neoformed bone (dark blue, in contact with a porous coated ceramic (+BMP-2) implanted in non osseous location (t= 3 months)

Luminescent nanoprobes for medical imaging

Nanostructured coatings on implants

Colloïdal apatite nanoparticles (20-100 nm), luminescent (rare earth element), associated with a targeting agent towards cancer cellules





Internalisation of uminescent nanoparticles l Within a human pancreatic cell



Application

Biomimetic apatite coating



HA- β TCP porous granules (Teknimed)

- with or without biomimetic apatite coating
- with or without adsorbed rh-BMP2





Application: Biomimetic apatite coating

Time (days)



H. Autefage et al. J. Biomed. Mater. Res. B (2009)



Application: Biomimetic apatite coating

In vivo study

Intramuscular implantation in sheep for 3 months



Osteoinduction potential of the nanocrystalline coating alone ... but WHY ?!! nanostructure ? nanoporosity ? composition ? Multiple parameters ?

... a question to investigate within NEWGEN network?

See Margarida Almeida talk (effect of [nanoHA] on ALP and BMP-2 expression)



Biomimetic apatite coating: Technology transfer

Coll. Teknimed S.A - CIRIMAT:

Ceraform®Revolution is a porous bioceramic with a nanocrystalline carbonated apatite coating enhancing surface reactivity leading to better bone reconstruction



H. Autefage et al. US patent n°12/487,101 (2009)



Conclusion

- Bone mineral is a very reactive complex structure
- The evolution of bone apatite with age corresponds to a physical-chemical maturation process
- Nanocrystalline apatites exhibit specific surface structure in aqueous media
- This hydrated surface layer is unstable
- Although the structure of the hydrated surface layer of apatites is unknown it offers ion exchange and protein adsorption capabilities
- The surface reactivity thus obtained is utilised by living organisms
- We can take advantage of this surface reactivity for making new bioactive bone substitute materials at low temperature

Tailoring calcium phosphate nanostructured materials for biomedical applications

Elizabete Costa (elisabete.costa@ua.pt) Margarida Almeida (margarida@ua.pt)



Centre for Research in Ceramics and Composite Materials

Department of Materials and Ceramic Engineering

UNIVERSITY OF AVEIRO



DRUG DELIVERY SYSTEMS DDS

PARTICLE REQUIREMENTS

Morphology accurately engineered

is a prerequisite for a predictable performance of **DDS**

Synthesis procedures able to prompt particle design with size and shape surgically tailored



OLOGY ISSUES

n



MP1301 NEWGEN WG1

...To exploit strategies based on chelant assisted synthesis for producing apatite nanoparticles targeted to biomedical devices (drug delivery systems + bone tissue growth)



MP1301 NEWGEN WG1

Chemical precipitation in presence of citrate





Citric acid $- x H^+ \longrightarrow$ Citrate ion



MP1301 NEWGEN WG1

Hydroxyapatite Synthesis Strategy

3 steps ...

Precursor solution Preparation

Solution heating

HAp particles



the lab

6

MP1301 NEWGEN WG1





MP1301 NEWGEN WG1





MP1301 NEWGEN WG1

Characterizing the precipitates ...

- Chemical Composition (ICP) and Crystallinity (X-ray Diffraction)
- Morphology (Electronic Microscopy)
- Specific surface area (BET)
- Electrical surface charge (Zeta potencial)
- Adsorbed species (FTIR)

Precipitation of HAp particles



MP1301 NEWGEN WG1

Nantes, 7th - 8th May 2014



[1] M.Martins et al, J.Colloid and Interface Sci, 2008;

[2] C.Santos et al, Microscopy and Microanalysis, 2009



MP1301 NEWGEN WG1





[Citrate]

pH=8,2

125nm

[1] M.Martins et al, J.Colloid and Interface Sci, 2008;

[2] C.Santos et al, Microscopy and Microanalysis, 2009



MP1301 NEWGEN WG1



55 m²/g



Cit:Ca=3:1

] M.Martins et al, J.Colloid and Interface Sci, 2008;



Hexagonal prisms



[2] C.Santos et al, Microscopy and Microanalysis, 2009





Nantes, 7th - 8th May 2014









WP1301 NE WG1 Nantes, 7th - 8th May 2014

Comparing HAp nanoparticle 55 m²/g SSA=170 m²/g







1] M.Martins et al, J.Colloid and Interface Sci, 2008;

HAp

[2] C.Santos et al, Microscopy and Microanalysis, 2009

nanoparticles : 180°C versus 37°C



MP1301 NEWGEN WG1

Nantes, 7th - 8th May 2014

Comparing HAP properties

HAp nanoparticles : 180°C versus 37°C







Santos, C. et al. Journal of the Royal Society Interface, 2012



MP1301 NEWGEN WG1

Interaction with cells

MG63 osteoblast-like cell response to HAp nanoparticles*



HAPnanoparticles concentration

Santos, C. et al. Journal of the Royal Society Interface, 2012



MP1301 NEWGEN WG1

Nantes, 7th - 8th May 2014



Low levels (not toxic) increased the expression of ALP and BMP-2

Santos, C. et al. Journal of the Royal Society Interface, 2012



Interaction with cells

MP1301 NEWGEN WG1

Nantes, 7th - 8th May 2014



NanoHap are readily internalized by MG63 cells. Toxic effects observed at high concentrations reflect the vesicular entrapment.

Santos, C. et al. Journal of the Royal Society Interface, 2012



MP1301 NEWGEN **WG1**

Concluding ...

HAP nucleation and growth phenomena are strongly dependent onprecursor solution chemistry { pH, [cit] } and temperature. Manipulatingthese variables allows controlling the shape, size and surface properties ofHAP nanoparticles (nanoHAp).

nanoHap are readily internalized by MG63 osteoblastic cells. Low [nanoHAp] do not affect cell viability/proliferation and enhance ALP and BMP-2 expression \rightarrow interesting profile for bone tissue applications.



MP1301 NEWGEN WG1

Acknowledgments

UNIVERSITY OF AVEIRO















Laboratório de Farmacologia e Biocompatibilidade Celular Permetentes é lineargatients

FCT Fundação para a Ciência e a Tecnologia MINISTÉRIO DA EDUCAÇÃO E CIÊNCIA